

THE SAFEST AND MOST ADVANCED

DISPOSABLE EVACUATION VALVES ON THE MARKET

Introducing the **NEWLY** designed **DOVE®** Disposable HVE and Saliva Ejector Backflow Prevention Valves. Alleviate cross contamination risks between patients and eliminate timely valve maintenance while saving time and money.

VIRUSES AND BACTERIA

survive for several days in the dark, moist environment of evacuation tubing, and disease transmission can occur when negative pressure causes the infectious material in the tubing to flow back into the next patent's mouth.*

ELIMINATE THE RISK

with **DOVE®** Disposable HVE and Saliva Ejector Backflow Prevention Valves. Accommodating most dental units, **DOVE®** one-way valves prevent backflow and eliminate cross contaminiation between patients.

*Recent research suggests the majority of dental offices DO NOT follow the mfg IFU for cleaning evacuation valves. Field swab bacteria surface testing concluded nearly 100% failure of providing a clean valve between patients.



OUR INNOVATIVE PRODUCT LINE

























Part Number | R2010-1-BL

HVE Assy Short Valve with Small Barb connection in blue color accommodates small barb connections. (DOVE® Model TP2000SM HVE Tailpiece Small Barb or similar tailpiece). Package Contains: Case Qty 600 (150/bag, 4 bags/case)

Part Number | R2010-1-PP

HVE Assy Short Valve with Small Barb connection in purple color accommodates small barb connections. (DOVE® Model TP2000SM HVE Tailpiece Small Barb or similar tailpiece). Package Contains: Case Qty 600 (150/bag, 4 bags/case)

Part Number | R2010-3-BL

HVE Assy Short Valve with Large Barb connection in blue color accommodates large barb connections. (DOVE® Model TP2000 Standard HVE Tailpiece or TP2000W HVE Tailpiece Wide or similar tailpiece).

Package Contains: Case Qty 600 (150/bag, 4 bags/case)

Part Number | R2010-3-PP

HVE Assy Short Valve with Large Barb connection in purple color accommodates large barb connections. (DOVE® Model TP2000 Standard HVE Tailpiece or TP2000W HVE Tailpiece Wide or similar tailpiece).

Package Contains: Case Qty 600 (150/bag, 4 bags/case)

Part Number | R2020-1-BL

HVE Assy Long Valve with Small Barb connection in blue color accommodates small barb connections. (DOVE® Model TP2000SM HVE Tailpiece Small Barb or similar tailpiece).

Package Contains: Case Qty 600 (150/bag, 4 bags/case)

Part Number | R2020-1-PP

HVE Assy Long Valve with Small Barb connection in purple color accommodates small barb connections. (DOVE® Model TP2000SM HVE Tailpiece Small Barb or similar tailpiece).

Package Contains: Case Qty 600 (150/bag, 4 bags/case)

Part Number | R2001B-BL

Saliva Ejector Backflow Valve in blue color. (Accomodates DOVE® Model TP2001 or similar tailpiece).

Package Contains: Case Qty 600 (150/bag, 4 bags/case)

Part Number | R2001B-PP

Saliva Ejector Backflow Valve in purple color. (Accomodates DOVE® Model TP2001 or similar

Package Contains: Case Qty 600 (150/bag, 4 bags/case)

Part Number | TP2000

The standard HVE Tailpiece securely connects to the HVE Tubing to accommodate DOVE® HVE Connection Valves (R2010-3-BL and R2010-3-PP styles).

Recommended For Dental Chair Hangers: Adec, Belmont, Forest, Midmark, Ritter, TPC.

Part Number | TP2000W

The wide HVE Tailpiece securely connects to the HVE Tubing to accommodate DOVE® HVE Connection Valves (R2010-3-BL and R2010-3-PP styles).

Recommended For Dental Chair Hangers: Dental EZ, Engle, Marus, Pelton and Royal/Promo.

Part Number | TP2000SM

The small barb HVE Tailpiece securely connects to the HVE Tubing to accommodate DOVE® HVE Connection Valves (R2020-1-BL, R2020-1-PP, R2010-1-BL, R2010-1-PP styles).

Recommended For: Most dental HVE hangers that accommodate the valve hanging by the neck or by the tailpiece.

Part Number | TP2001

The Saliva Ejector Tailpiece securely connects to the SE Tubing to accommodate DOVE® SE Connection Valves (R2001B-BL and R2001B-PP styles).